





Aesculap Power Systems

microspeed® uni

-
- (GB)** **Instructions for use**
microspeed® uni motor system
 - (USA)** **Instructions for use**
microspeed® uni motor system
 - (D)** **Gebrauchsanweisung**
microspeed® uni-Motorensystem
 - (F)** **Mode d'emploi**
➤ voir www.aesculap-extra.net
 - (E)** **Instrucciones de manejo**
➤ ver www.aesculap-extra.net
 - (I)** **Istruzioni per l'uso**
➤ vedere www.aesculap-extra.net
 - (P)** **Instruções de utilização**
➤ ver www.aesculap-extra.net
 - (NL)** **Gebruiksaanwijzing**
➤ zie www.aesculap-extra.net
 - (S)** **Bruksanvisning**
➤ se www.aesculap-extra.net
 - (FIN)** **Käyttöohjeet**
➤ ks. www.aesculap-extra.net
 - (RUS)** **Инструкция по применению**
➤ смотри www.aesculap-extra.net

microspeed® uni motor system

Legend

- Captions first foldout page: Control unit
- 1 microspeed® uni control unit with cooling unit GD670
 - 2 Display and touch control field
 - 3 Release buttons
 - 4 Coolant pump
 - 5 Motor connection sockets
 - 6 Power OFF switch
 - 7 Power ON switch
 - 8 Power ON indicator
 - 9 Bottle holder
 - 10 Bottle holder mount
 - 11 Ventilation grate
 - 12 Power socket
 - 13 Fuse holder
 - 14 Connection socket for foot control
 - 15 Stacking cones
 - 16 Locking screw
 - 17 Equipotentialization connector
 - 39 Interfaces to be used by the manufacturer only
- Captions second foldout page: Foot controls, motor cables, motors
- 18 Foot control function button
 - 19 Right pedal
 - 20 Left pedal
 - 21 microspeed® uni foot control, two pedals, GD671
 - 22 Motor direction switch
 - 23 Pedal
 - 24 microspeed® uni foot control, one pedal, GD668
 - 25 Motor release switch
 - 26 Hand control function button
 - 27 microspeed® uni motor cable with hand control GD673
 - 28 Release key
 - 29 Lever
 - 30 Motor release switch
 - 31 microspeed® uni motor cable for foot control GD672
 - 32 Motors with motor-side plug connection
 - 33 microspeed® uni mini pistol handpiece GD684
 - 34 Upper trigger
 - 35 Lower trigger
 - 36 Suction control lever
 - 37 microspeed® uni shaver handpiece GD682
 - 38 Tube olive

Contents

- 1. Safe handling 3
- 2. Product description 4
 - 2.1 System components 4
 - 2.2 Components necessary for using the product 4
 - 2.3 Intended use 4
 - 2.4 Mode of operation 4
 - 2.4.1 Control unit 4
 - 2.4.2 Control and display interface 4
 - 2.4.3 Control fields on the display 5
 - 2.4.4 Motor control fields 5
 - 2.4.5 Pump control field 6
 - 2.4.6 Instruction menu 6
 - 2.4.7 Motor recognition 6
 - 2.4.8 Automatic saving of device settings 6
 - 2.4.9 Overload cut-out 7
 - 2.4.10 Coolant pump 7
 - 2.4.11 Motor cable 7
 - 2.4.12 Motors 9
 - 2.4.13 Foot control devices 11
 - 2.5 Symbols 12
 - 2.5.1 Symbols on the control unit 12
 - 2.5.2 Symbols on foot control, motor cable and motor 12
 - 2.5.3 Symbols in display 13
- 3. Preparation and setup 15
 - 3.1 Preparation 15
 - 3.1.1 Mounting the control unit on the mobile stand 15
 - 3.1.2 Dismounting the control unit from the mobile stand 15
 - 3.1.3 Mains power connection 15
- 4. Working with the microspeed® uni motor system 16
 - 4.1 System set-up 16
 - 4.1.1 Connecting the accessories 16
 - 4.1.2 Connecting the foot control to the control unit 16
 - 4.1.3 Connecting the tubing set 17
 - 4.1.4 Connecting the motor cable to the control unit 17
 - 4.1.5 Connecting a motor (GD674, GD676, GD677, GD678 and GD679) to the motor cable 18
 - 4.1.6 Switching on the control unit/Automatic self-test 18
 - 4.1.7 Connecting handpieces/attachments/blades 18
 - 4.2 Function checks 22
 - 4.3 Safe operation 23
 - 4.3.1 Setting up and operating the motors GD674, GD676, GD677, GD678 and GD679 23
 - 4.3.2 Setting up and operating the microspeed® uni shaver handpiece GD682 26
 - 4.3.3 Setting up and operating the microspeed® uni mini pistol handpiece GD684 28
 - 4.4 Settings in the Instruction menu 32
 - 4.4.1 Submenu information on the error conditions shown on the display 33
 - 4.4.2 Submenu Device-specific settings 33
 - 4.4.3 Submenu Motor (type)-specific settings 35

5.	Processing	38
5.1	Processing the microspeed® uni control unit with cooling unit GD670 and microspeed® uni foot control devices GD668/GD671	38
5.1.1	Cleaning/Disinfection	38
5.1.2	Manual cleaning/disinfecting	39
5.1.3	Control, care and inspection	39
5.1.4	Storage	39
5.2	Processing the microspeed® uni motors GD674, GD676, GD677, GD678, GD679 and motor cables GD672 and GD673	39
5.2.1	Preparations	39
5.2.2	Cleaning/Disinfection	39
5.2.3	Mechanical cleaning/disinfecting	40
5.2.4	Manual cleaning/disinfecting	40
5.2.5	Control, care and inspection	40
5.2.6	Packaging	40
5.2.7	Sterilization method and parameters	41
5.2.8	Sterilization for the US market	41
5.2.9	Storage	41
5.2.10	Maintenance	41
5.3	Processing the microspeed® uni shaver handpiece GD682	41
5.3.1	Preparations	41
5.3.2	Mechanical cleaning/disinfecting	41
5.3.3	Manual cleaning/disinfecting	42
5.3.4	Control, care and inspection	42
5.3.5	Packaging	42
5.3.6	Sterilization method and parameters	42
5.3.7	Sterilization for the US market	42
5.3.8	Storage	42
5.3.9	Maintenance	43
5.4	Processing the microspeed® uni mini pistol handpiece GD684	43
5.4.1	Preparations	43
5.4.2	Cleaning/Disinfection	43
5.4.3	Mechanical cleaning/disinfecting	43
5.4.4	Manual cleaning/disinfecting	44
5.4.5	Control, care and inspection	44
5.4.6	Packaging	44
5.4.7	Sterilization method and parameters	44
5.4.8	Sterilization for the US market	44
5.4.9	Storage	44
5.4.10	Maintenance	45
6.	Maintenance	45
7.	Troubleshooting list	46
7.1	Malfunctions with error numbers	46
7.2	Other malfunctions	48
7.3	Fuse change	49
8.	Technical Service	49
9.	Accessories/Spare parts	49
9.1	Accessories	49
9.1.1	microspeed® uni motors, motor cables and foot control devices	49
9.1.2	Cooling unit	49
9.1.3	Other components	49
9.1.4	Power cord	50
9.1.5	Maintenance, processing and service	50

9.2	Spare parts	50
10.	Technical specifications	50
10.1	Classification acc. to directive 93/42/EEC	50
10.2	microspeed® uni control unit with cooling unit GD670	51
10.2.1	Factory settings	51
10.3	Low-speed motors	51
10.3.1	Low-speed motors with micro-Line coupling	51
10.3.2	Low-speed motor with mini-Line coupling	52
10.4	High-speed motors	52
10.4.1	High-speed motors with Hi-Line coupling	52
10.5	microspeed® uni shaver handpiece GD682	52
10.6	microspeed® uni mini pistol handpiece GD684	53
11.	Disposal	53
12.	Distributor in the US/Contact in Canada for product information and complaints	53

1. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Risk of death by electric shock!

- The user must not open the control unit under any circumstances.
- Connect product only to mains with protective circuit.

Note

For handling and operating the products that are used with the microspeed® uni motor system, always observe the instructions for use of the respective products. The present instructions for use only describe the handling and operation of the "microspeed® uni motor system".

- Prior to use, check for proper condition and functioning of the product.
- Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.
- To avoid damaging the product through improper setup or operation and mitigating the warranty and liability on the part of the product's manufacturer:
 - Use the product only according to these instructions for use.
 - Follow the safety instructions and maintenance advisories.
 - Only combine Aesculap products with each other.
- Ensure that the products and accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Keep the instructions accessible for the O.R. personnel.
- Always adhere to applicable standards.
- Verify that the electrical installations in the room in which the equipment is to be used conform to IEC standards.
- Unplug the device by pulling on the plug, and never on the power cord.
- Do not use the control unit and motors in explosion hazard areas.
- Prior to use, carry out sterile processing of motor(s), motor cable, tubing, handpiece(s) and tools.
- Operate the handpieces only at the recommended motor speeds.

microspeed® uni motor system

2. Product description

2.1 System components

Designation	Art. no.
microspeed® uni control unit with cooling unit	GD670
Bottle holder	GD412804
Instructions for use	TA022089
Notes on Electromagnetic Compatibility (EMC)	TA022130

2.2 Components necessary for using the product

- microspeed® uni control unit with cooling unit GD670
- Power cord
- Motor
- Motor cable
(when using a motor with a motor-side plug connection)
- Foot control
(when using a motor cable without hand control, or the shaver handpiece)
- Bottle holder
(when using the coolant pump)
- Tubing set
(when using the coolant pump)
- Surgical suction device
(when using the shaver handpiece)

2.3 Intended use

The microspeed® uni motor system is a universal system with respect to its range of applications and to the system composition and its variety of uses.

The microspeed® uni motor system is used in orthopedics/traumatology, neurosurgery, spine surgery, hand/foot surgery, ENT/OMF surgery, plastic surgery and arthroscopy.

The microspeed® uni motor system can be operated as complete from the sterile area, either through foot control or through hand control.

2.4 Mode of operation

The microspeed® uni motor system comprises the following components:

- Control unit with coolant pump, for universal application
- Two motor cables (with or without hand control)
- Various motors with different coupling types
- Pistol handpiece
- Shaver handpiece
- Two foot control devices

2.4.1 Control unit

Control unit 1 has been designed for a mains voltage rating of 100 V to 240 V/50 Hz to 60 Hz. The mains voltage is transformed so that the motors run at safety low voltage.

Control unit 1 is equipped with two motor connection sockets 5 for connecting two (different) motors, and with one foot control connection socket 14 (on rear panel).

The following devices can be connected to control unit 1 at the same time:

- Two motors and
- One foot control (GD668 or GD671).

Note

Do not connect two shaver handpieces GD682 at the same time.

Only one motor can be activated at any one time.

Motors with motor-side plug connection (GD674, GD676, GD677, GD678 and GD679) can be operated either through hand control (motor cable GD673) or foot control GD668 or GD671 (with motor cable GD672).

Note

When two motors without hand control are connected to the control unit (motor cable GD672), a motor can only be activated if the motor release switch of the motor that is to be activated is in the "ON" position while the other motor is locked (motor release switch in "OFF" position).

2.4.2 Control and display interface

Display 2 shows the present device status (operating and error conditions) at all times. The display is divided into motor control field(s) and a pump control field. The device parameter settings can be changed by pressing the control keys.

For information on all the symbols used in the motor and pump control fields, see chapter Symbols in display.

Control key "Instruction menu" calls up the Instruction menu.

2.4.3 Control fields on the display


Fig. 1

Legend

- A Motor control field(s)
- B Pump control field
- C Instruction menu

The display is divided into three different control fields, which are described in the following chapters.

2.4.4 Motor control fields


Fig. 2

Legend

- A Status display and control key for releasing (symbol highlighted by green background)/locking the pump
- B Activation mode indicator: hand control or foot control
- C Present motor speed range upper limit setting
- D Qualitative display of the actual motor speed
- E Control keys for increasing or reducing the motor speed range upper limit
- F Motor direction indicator
- G Operating mode selected for the connected motor type, e.g. saw program for motors GD678 and GD679
- H Connected motor type
- I Maximum selectable motor speed range upper limit for the connected motor type or with the selected operating mode, respectively

Each motor connected to control unit 1 has a motor control field assigned to it.

Two motor control fields are opened in display 2 only if two motors 32 are actually connected to control unit 1.

Only one motor control field will be opened if only one motor is connected.

The control fields are visually identified as belonging to the respective motor connection sockets 5.

microspeed® uni motor system

A motor control field is displayed as an active control field when the motor assigned to the motor control field:

- is operated
- was the last motor operated.

The symbol of the connected motor type is highlighted by a green background whenever a motor control field is released for operation.

2.4.5 Pump control field


Fig. 3

Legend

A Status indicator symbol of coolant pump 4. The symbol is highlighted by a green background while the pump is active. The pump can only be activated if it has been released for operation and if the motor is active.

B Control keys for increasing or reducing the delivery rate

C Display of selected dosage

The preset dose level **C** applies to the motor whose settings are displayed in the (white background) motor control field.

2.4.6 Instruction menu

Note

The Instruction menu can only be called up at the control unit.

This control key is pressed to call up the Instruction menu. In the Instruction menu, information about the error conditions output to the display and device settings such as the user language and the volume of the acoustic alarm/notification signal can be changed.

Motor-type specific adjustments and settings can be carried out here, too:

- Slow-down and speed-up ramps for each motor type
- Oscillating or tapping mode (for pistol handpiece GD684)
- Oscillation frequency and speed-up ramp in oscillating mode (in Shaving mode of operation of shaver handpiece GD682)
- Reset: to restore the factory settings

2.4.7 Motor recognition

Control unit 1 recognizes the motor type connected to one of the motor connection sockets 5 of control unit 1 and displays the type in the respective motor control field of display 2. The most recent settings for this motor type at the same motor connection socket are shown on display 2.

The settings can be changed within the motor-specific setting ranges, either on the display or through hand or foot control (in the Device settings menu), see chapter Technical specifications.

The motor control field of the most recently activated motor is always brightly highlighted, whereas the motor control field of the motor that is not the one that was most recently activated or locked is displayed with a gray background. When only one motor is connected to control unit 1, the corresponding control field is highlighted, if this motor is not locked.

Saw program (only for low-speed motors GD678 and GD679)

In this program the maximum motor speed is limited to 16,000 1/min. The settings most recently used with the saw program (motor speed range upper limit, dose level for coolant pump 4) are recalled automatically.

2.4.8 Automatic saving of device settings

The most recent settings for this motor type at this motor connection socket (speed range upper limit, delivery rate, etc.), are called up automatically as soon as the same motor type is recognized again as connected at the same motor connection socket of control unit 1 after a motor change.

Note

After each device start-up both the coolant pump and motor type-specific special functions such as the saw program (for GD678 and GD679) or the oscillating or tapping mode (for GD684) are always inactive (even if they had been released or activated before). All parameters are reset when the operating key "Factory setting" in the Instruction menu is pressed.

2.4.9 Overload cut-out

To protect the connected motors against thermal destruction, the system features a separate motor temperature monitoring device. For this purpose, a temperature module is integrated in the device.

- 1st Level:
As soon as the temperature module indicates a motor temperature above 75 °C and below 90 °C, an acoustic warning is sounded and the message "Temp. Motor" flashes on the display. Any further overload causes control unit 1 to switch to level 2.
- 2nd Level:
As soon as a motor temperature in excess of 90 °C is detected, the motor is switched off and the message "Error 7 – Motor overheated" is shown on display 2.

After a cooling-down period the motor can be activated again.

Alternatively, the temperature module is reset, and work can be continued as soon as the overheated motor has been replaced. We recommend keeping a replacement motor on standby.

2.4.10 Coolant pump

Note

The coolant pump only works when a motor is running!

At power-on of the control unit, coolant pump 4 is always deactivated (even if it had been released for operation at an earlier time).

The microspeed® uni control unit 1 is fitted with a coolant pump 4. This coolant pump serves exclusively for cooling or irrigating the tools. The delivery rate at the highest dose level is approx. 80 ml per minute. The coolant pump can be activated either in the respective motor control field or through the appropriate hand or foot control device. The delivery rate can be selected via the pump control field or through the Device settings menu, see chapter Submenu Device-specific settings.

Exception: When using pistol handpiece GD684, the pump can be activated, and the delivery rate can be changed, only via the appropriate control fields on the control unit.

Any one of 21 dosage settings can be selected. At dosage settings 1 to 5, coolant pump 4 operates in intermittent mode. The pump can be switched on and off, and the delivery rate can be changed, even while the motor is active. Once the motor was active for more than two seconds, a special drive control prevents further coolant from dripping out.

2.4.11 Motor cable

Risk of injury due to uncontrolled movement of skin abrading handpiece GB280 in conjunction with a motor cable with hand control GD673!

- **Use skin abrading handpiece GB280 only with motor cable GD672 and foot control.**

Note

A motor cable is needed to connect motors (with motor-side plug connection) GD674, GD676, GD677, GD678 and GD679 to the control unit.

The microspeed® uni motor system comprises one motor cable for foot control (GD672) and one motor cable with hand control (GD673). Both motor cables can be used in connection with all motors with a motor-side plug connection. Either of the two motor cables can be connected at any of the two motor connection sockets 5.

Notes

The choice of motor cable depends on the application and which handpiece is to be used.

For using saws GB128R, GB129, GB130R and motor cable with hand control GD673, we recommend guiding the saw with both hands for optimal handling. If single-handed saw guidance is required, use motor cable GD672 with a foot control device.

microspeed® uni motor system

microspeed® uni motor cable for foot control GD672


Fig. 4

Legend

30 Motor release switch

Note

This motor cable can only be used in connection with a foot control.

The motor coupled to this motor cable can only be activated if motor release switch **30** is in the "ON" position. When two motor cables GD672 are used at the same time, a motor can be activated only if it has been released (motor release switch **30** in position "ON") while the other motor is deactivated (motor release switch **30** in position "OFF").

microspeed® uni motor cable with hand control GD673


Fig. 5

Legend

25 Motor release switch:

- Position "ON": Motor released for operation
- Position "OFF": Motor locked against operation

26 Hand control function button:

- Short actuation: Switching on/off the coolant pump
- Pressing for more than three seconds (with the motor inactive): Change to Device settings menu, siehe Kapitel Changing the settings in the Device settings menu through hand control

29 Lever:

- Actuator for motor speed setting

Note

This motor cable can not be used in connection with a foot control.

The motor coupled to this motor cable can only be activated if the motor release switch is in its "ON" position.

2.4.12 Motors

Motors with motor-side plug connection

Note

A motor cable GD672 or GD673 is required for connecting motors with motor-side plug connection to the control unit!

Depending on the field of application, there are five different motor types available, with the three Aesculap standard coupling systems Hi-Line, micro-Line (Intra coupling system to DIN 13940/ISO 3964) and mini-Line. Each motor can be operated in right or left (clockwise or counterclockwise) rotation mode. When a motor is running in left-rotation mode, control unit 1 issues an acoustic signal.

For the performance data and parameters to be set for the individual motor types, see chapter Technical specifications.

Motor types GD678 and GD679 offer the option to select a saw program. In this program the maximum motor speed is limited to 16,000 1/min.

Motors with fixed motor cable GD682 and GD684

Note

The motors are directly connected to the control unit. An additional motor cable is not required!

For the performance data and parameters to be set for the individual motor types, see chapter Technical specifications.

microspeed® uni shaver handpiece GD682

Fig. 6

Legend

36 Suction control lever:

Lever for continuous regulation of the suction power

38 Tube olive:

Connection to the surgical suction device via a suction tube

Notes

Always observe the instructions for use of blades TA011579.

Only use Aesculap blades.

Shaver handpiece GD682 can only be used with foot control GD671 and a surgical suction device (e.g. GF200 or GF210).

Shaver handpiece GD682 can be used if:

- a motor with motor cable GD673 is connected at the second motor connection socket, or
- a pistol handpiece GD684 is connected at the second motor connection socket, or
- a motor with motor cable GD672 is connected at the second motor connection socket and the motor release switch of that cable is in its "OFF" position, or
- only one shaver handpiece GD682 is connected to control unit 1.

There are two types of blades available for shaver handpiece GD682: reamer and shaver blades.

Control unit 1 recognizes whether the coupled blade is a reamer or a shaver blade. The respective mode of operation (reaming or shaving) is indicated on the display.

An acoustic signal is issued whenever a blade of a different type is inserted. At the same time the respective symbol is displayed in the appropriate control field.

With a reamer blade inserted, left-rotation and right-rotation modes are available, but not the oscillating mode. The maximum speed is limited to 8,000 1/min. The symbol for the "Reaming" operating mode is displayed in the control field.

With a shaver blade inserted, right-rotation and left-rotation operating modes are available. The maximum speed is limited to 5,000 1/min. The symbol for the "Shaving" operating mode is displayed in the control field.

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- 34 Upper trigger
- 35 Lower trigger

Note

The pistol handpiece can be locked against inadvertent actuation by gently pressing and turning the upper trigger by 90°.

The pistol handpiece cannot be used with foot control.

Pistol handpiece GD684 is fitted with the Aesculap mini-Line coupling system.

A protective sleeve is supplied with the handpiece. This sleeve serves to protect against injury and damage when inserting Kirschner wires.

There are two operating modes of pistol handpiece GD684:

- Standard operation
- Oscillating or tapping mode

If control unit 1 recognizes a pistol handpiece GD684 connected, the pistol handpiece will operate in the Standard mode (factory setting). The device parameter settings (speed range upper limit, delivery rate, pump status) can be changed only via the touch control field of control unit 1. Oscillating or tapping mode can be selected via the Instruction menu. The symbol for "Oscillating or tapping mode" is displayed in the respective motor control field. Oscillating or tapping mode is activated through the upper trigger. An acoustic notification signal is sounded and "Oscillating or tapping mode active" flashes in the respective motor control field. In this mode of operation, the speed range upper limit cannot be changed via the touch control field on control unit 1.

The oscillating or tapping mode can be deactivated again by pressing the upper trigger. After each power-up of the control unit and connection of the pistol handpiece, oscillating or tapping mode is always deactivated and needs to be activated anew via the pistol handpiece. The preselected setting in the Instruction menu is saved and remains in memory even through subsequent system restarts.

Standard operation:

Right and left rotation mode (when in left-rotation mode, control unit 1 issues an acoustic signal).

Motor direction switch/trigger for left rotation 34, see Fig. 7:

- Pressing lower trigger 35 while holding down upper trigger 34 makes the pistol handpiece operate in left-rotation mode.
- The motor direction is changed by pressing trigger 34 with trigger 35 already actuated.

Motor speed control trigger 35:

- Actuator for motor speed setting (motor speed proportional to actuator travel).
- If only trigger 35 is pressed, the pistol operates in right rotation mode.

Oscillating or tapping mode:

Activating or deactivating the operating mode, see Fig. 7:

- This mode of operation is activated (and deactivated) by pressing and holding upper trigger 34 for three seconds.

Working in "Oscillating or tapping mode":

- The pistol handpiece is operated in tapping mode by pressing lower trigger 35 within the first half of the total trigger travel. The motor speed is proportional to the actual trigger travel.

Tapping mode means: Clockwise rotation by a certain angle (forward step) followed by counterclockwise rotation (backward step) by a smaller angle.

- The pistol handpiece is operated in oscillating mode by pressing lower trigger 35 within the second half of the total trigger travel. The motor speed is proportional to the actual trigger travel.

Oscillating mode means: Clockwise rotation by a certain angle followed by counterclockwise rotation by the same angle.

The pistol handpiece is operated in left-rotation mode by pressing lower trigger 35 at the same time as upper trigger 34.

2.4.13 Foot control devices

Notes

The *microspeed® uni* motor system includes two different foot controls (GD668 and GD671).

Only one foot control can be connected to control unit GD670 at any time.

microspeed® uni foot control, one pedal, GD668

Fig. 8

Legend

- 18** Foot control function button:
Short actuation: Switching on/off the coolant pump
Actuation for more than three seconds (with the motor inactive):
Changing to Device settings menu, see chapter Changing the settings in the Device settings menu through the foot control
- 22** Motor direction switch:
Selection of right or left rotation
- 23** Pedal:
Actuator for motor speed setting

microspeed® uni foot control, two pedals, GD671


Fig. 9


Legend

- 18** Foot control function button:
Short actuation: Switching on/off the coolant pump
Actuation for more than three seconds (with the motor inactive):
Change to Device settings menu, see chapter Changing the settings in the Device settings menu through the foot control.
- 19** Right pedal:
Actuator for motor speed setting in right rotation
- 20** Left pedal:
Actuator for motor speed setting in left rotation


microspeed® uni motor system

2.5 Symbols

2.5.1 Symbols on the control unit

	Adhere to instructions for use
	Type BF application component
	Input from irrigation bottle (coolant pump)
	Outlet to handpiece (coolant pump)
	Pumping direction (coolant pump)
	Equipotentialization connector
	Fuse
	Alternating current
	Foot control
	Marking of electric and electronic devices according to directive 2002/96/EC (WEEE), see chapter Disposal.

2.5.2 Symbols on foot control, motor cable and motor

	Foot control function button (Foot control GD668, GD671, motor cable GD673)
	Motor direction switch Selection of left or right rotation (Foot control GD668)
	Pedal for left/right rotation (according to the motor direction set through the motor direction switch, foot control GD668)
	Pedal for right rotation (Foot control GD671)
	Pedal for left rotation (Foot control GD671)
	Oscillating operation when both pedals are pressed (only allowed in "Shaving" mode in conjunction with shaver handpiece GD682 , foot control GD671)
	Class AP anesthetic test (EN 60601-1) (Foot controls GD668 and GD671)
IPX8	Enclosure protection type to DIN EN 60529 (Foot controls GD668 and GD671)
	Release (Motors GD674, GD682 and GD684)
	Motor speed control (Motor GD684)
	Right/left rotation switch (Motor GD684)
	Oscillating or tapping mode: (Motor GD684)
	Protection against inadvertent activation (Motor GD684)

2.5.3 Symbols in display

Display symbols in the motor control fields

	(Set) motor speed range upper limit (example: 20,000 1/min)
	Motor speed in 1/min
	Qualitative display of the actual motor speed (for shaver handpiece GD682: actual tool speed) Symbol displayed on gray background in motor control field if the respective motor is not activated
	Qualitative display of the actual motor speed (for shaver handpiece GD682: actual tool speed); active while motor is operated
	Motor direction indicator for right rotation when the motor is inactive or when right rotation was not selected through foot control GD668; in an inactive motor control field, this symbol is displayed on a gray background
	Motor direction indicator for right rotation active when the motor is operated or in standby, provided this motor direction was selected and foot control GD668 is connected to the control unit
	Motor direction indicator for left rotation when the motor is inactive or when left rotation was not selected through foot control GD668; in an inactive motor control field, this symbol is displayed on a gray background
	Motor direction indicator for left rotation active when the motor is operated or in standby, provided this motor direction was selected and foot control GD668 is connected to the control unit
	Hand control active if connected hand control is recognized
	Foot control active if connected foot control is recognized and at least one motor is not operated through a hand control device
	No response when foot control is pressed down active if neither a foot control nor a hand control device is recognized as connected
	GD676 microspeed® uni Hi 150 motor (high-speed motor; 150 W; coupling system: Hi-Line (XS) connected) If the motor can be activated (is released), the symbol is highlighted by a green background

	GD677 microspeed® uni Hi 100 motor (high-speed motor; 100 W; coupling system: Hi-Line (XS) connected) If the motor can be activated (is released), the symbol is highlighted by a green background
	GD678 microspeed® uni micro 150 motor (low-speed motor; 150 W; coupling system: micro-Line) If the motor can be activated (is released), the symbol is highlighted by a green background
	GD679 microspeed® uni micro 100 motor (low-speed motor; 100 W; coupling system: micro-Line) If the motor can be activated (is released), the symbol is highlighted by a green background
	GD674 microspeed® uni mini 100 motor (low-speed motor; 100 W; coupling system: mini-Line) connected If the motor can be activated (is released), the symbol is highlighted by a green background
	GD682 microspeed® uni shaver handpiece connected If the motor can be activated (is released), the symbol is highlighted by a green background
	GD684 microspeed® uni mini pistol handpiece (coupling system: mini-Line) connected If the motor can be activated (is released), the symbol is highlighted by a green background
	Operating mode "Reaming" (only in connection with shaver handpiece GD682, if no shaver tip is recognized)
	Operating mode "Shaving" (only in connection with shaver handpiece GD682, if a shaver tip is recognized)
	Oscillating or tapping mode (only in connection with pistol handpiece GD684) Displayed when oscillating or tapping mode was released in the instruction menu, but not activated through the pistol handpiece

microspeed® uni motor system


	Oscillating or tapping mode active: (only in connection with pistol handpiece GD684) Displayed when oscillating or tapping mode was released in the Instruction menu and is activated through the pistol handpiece (by holding down the upper trigger for three seconds)
---	--

	"Pump ON" indicator and control key for activating the pump Active when, in connection with the selected motor, the pump is switched on via the function button (hand or foot control) or by actuating this field. The symbol is displayed in the respective motor control field.
---	--

Display/Control key symbols in the motor control fields

	Increasing the motor speed range upper limit
	Reducing the motor speed range upper limit
	Status indicator/control key of the saw program This symbol indicates that the saw program, which can only be selected in connection with motors GD678 and GD679, is inactive
	Status indicator/control key of the saw program This symbol indicates that the saw program, which can only be selected in connection with motors GD678 and GD679, is active
	Motor direction indicator when left rotation was not selected, and control key for selecting left rotation The control key symbol is only displayed, and a motor direction can only be selected if hand control GD673 is connected
	Motor direction indicator when left rotation was selected, and control key for selecting left rotation The control key symbol is only displayed, and a motor direction can only be selected if hand control GD673 is connected
	Motor direction indicator when right rotation was not selected, and control key for selecting right rotation The control key symbol is only displayed, and a motor direction can only be selected if hand control GD673 is connected
	Motor direction indicator when right rotation was selected, and control key for selecting right rotation The control key symbol is only displayed, and a motor direction can only be selected if hand control GD673 is connected
	"Pump OFF" indicator and control key for activating the pump The symbol is displayed in the respective motor control field

Display symbols in the pump control field

	"Pump inactive" indicator Symbol as displayed with the motor inactive, i.e. when the pump has not been activated in the pump control field
	"Pump active" indicator is active when the pump released in connection with the selected and activated motor is activated
	Display of selected dosage is always active

Control key symbols in the pump control field

	Increasing the delivery rate
	Reducing the delivery rate

Control key symbols within the display

	Selecting the secondary control level Instruction menu, can only be activated when no motor is active
---	---

3. Preparation and setup

Non-compliance with these rules will result in complete exclusion of any responsibility and liability on the part of Aesculap.

- Observe the following regulations when setting up and operating the product:
 - the national regulations for installation and operation
 - national regulations on fire and explosion protection.
 - operating advisories acc. to IEC-/VDE regulations.

The following advice is given on the basis of the presently applicable IEC-/VDE regulations:

Fire and explosion hazard!

- Do not use the control unit and motors in explosion hazard areas.

The foot control circuit is ignition-safe and approved for operation in medical environments acc. to IEC 601/VDE 0750.

3.1 Preparation

Note

The control unit can be mounted on an Aesculap mobile stand for *microspeed® uni (adjustable height) GA411* or on an Aesculap mobile stand for surgical devices *GD416M*.

Other Aesculap devices must not be stacked on top of the control unit.

3.1.1 Mounting the control unit on the mobile stand

Note

Follow the instructions for use of the mobile stand.


Fig. 10

Legend

- A Device support
- 15 Stacking cones
- 16 Locking screw

Only one device may be mounted on the mobile stand at any time, always taking into consideration the stability of the mobile stand.

- Remove the caps from stacking cones **15**.
- Remove the device feet, using a screwdriver.
- Move stacking cones **15** to the correct mounting positions by turning locking screw **16** counterclockwise as far as it will go, using a screwdriver.
- Position the device on the mobile stand in such a way that front edge sits flush with device support **A** and stacking cones **15** sit in the corresponding holes in device support **A**.
- Secure the device by turning locking screw **16** clockwise as far as it will go.
- Gently lift the device to check if it is securely attached to device support **A**.
The device is now mounted on the mobile stand.

3.1.2 Dismounting the control unit from the mobile stand

- Use a screwdriver to turn locking screw **16** counterclockwise as far as it will go.
- Lift off the device.
- Reattach the device feet either by hand or with a screwdriver.
The device is now dismounted from the mobile stand.

3.1.3 Mains power connection

The voltage of the mains power must correspond to the voltage indicated on the type plate of the device.

- Plug in the power cord at power cord socket **12**.
- Plug in the mains plug at the building mains power socket.

microspeed® uni motor system

4. Working with the microspeed® uni motor system

Risk of infection and contamination!

- Prior to first operation, carry out sterile processing of the motors, motor cable, tubing, handpiece and tools.
- Make certain that the packaging of sterile products is undamaged.

Risk of injury caused by inadvertent activation of the motor!

- Secure motors that are not actively operated against inadvertent activation through foot or hand control.

4.1 System set-up

4.1.1 Connecting the accessories

Note

Follow the instructions for use of individual accessories!

Note

micro-Line drill pistol (GB200) and arthroscopy handpiece (GB801) must not be used with this product.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically designed for the intended application and do not compromise the performance and safety characteristics of the products.

Depending on the motor type, all handpieces and tools of this range can be coupled through the three Aesculap standard coupling systems Hi-Line (or Hi-Line XS), micro-Line (handpieces for drilling, reaming, sawing and shaping, with Intra coupling system acc. to DIN 13940/ISO 3964) and mini-Line.

For other accessories, see chapter Accessories/Spare parts and brochure no. O22711, Aesculap Power Systems.

4.1.2 Connecting the foot control to the control unit

Note

Only necessary when using motor cable for foot control GD672 or shaver handpiece GD682.

The foot control can not be used in combination with motor cable with hand control GD673 and/or pistol handpiece GD684.


Fig. 11

Legend

- A** Foot control connector
- B** Locking pins
- C** Locking clips
- D** Arrows

14 Connection socket for foot control

- To connect the foot control connector **A** to control unit **1**: Fully plug in the foot control connector **A** at the connection socket for foot control **14** (in rear panel of the unit). Locking clips **C** audibly engage in locking pins **B**, see Fig. 11.
- To disconnect the foot control device from control unit **1**: Simultaneously compress locking clips **C** (see arrows **D**) and unplug connector **A**.

4.1.3 Connecting the tubing set

Note

The tubing set is sterile! Sterile separation is carried out at the tubing.


Fig. 12

Legend

- A Arrows
- B Arrows
- 3 Release buttons
- 4 Coolant pump
- To open coolant pump 4: Simultaneously press release buttons 3.
- Insert tubing: Insert the pump tube with the tube nipple into the recess (arrows B).
- Close coolant pump 4: Push together the pump housing until it clicks into position (arrows A). When doing this, ensure that the pump enclosure is completely closed, both at the inlet (from the irrigation bottle) and at the outlet (to the handpiece).

- Insert bottle holder 9 for the sterile fluid bottle into bottle holder mount 10.
- Insert the spike into the sterile fluid bottle.
- Mount the sterile fluid bottle at bottle holder 9.

4.1.4 Connecting the motor cable to the control unit

Note

The motor cable and the motor are sterile! Sterile separation is carried out at the motor cable.


Fig. 13

Legend

- A Mark (motor connection socket)
- B Plug
- C Mark (plug)
- D Arrows
- 5 Motor connection socket

Damage to, or destruction of, the motor cable caused by pulling the cable!

CAUTION

➤ Always pull at the plug, never at the cable.

- To connect plug B: Align plug B in such a way that mark C on the plug is aligned with mark A on motor connection socket 5.
- Fully insert plug B in the motor connection socket.
- To disconnect plug B: Hold plug B in the marked places (arrows D) and unplug it from motor connection socket 5.

microspeed® uni motor system

4.1.5 Connecting a motor (GD674, GD676, GD677, GD678 and GD679) to the motor cable


Fig. 14

Legend

- A Mark on motor cable
- B Mark on motor
- C Arrows
- 28 Release key (x2)
- 32 Motor

- To couple motor 32: Install motor 32 on the motor cable, making certain that mark B on the motor is aligned with mark A on the motor cable. Motor 32 will audibly engage in its proper position.
- To uncouple motor 32: Compress both release keys 28 at the points indicated by arrows C. (The motor will be pushed away slightly from the motor cable.) This completes the uncoupling of motor 32.

4.1.6 Switching on the control unit/Automatic self-test

- Switch on control unit 1 at mains power ON switch 7. Power ON indicator 8 will be illuminated. Control unit 1 always performs a self-test when the unit is switched on. During this test the microspeed® uni start picture is displayed. If any malfunction is detected, an error message will be displayed, see chapter Malfunctions with error numbers.

4.1.7 Connecting handpieces/attachments/blades

WARNING Risk of injury due to uncontrolled movement of skin abrading handpiece GB280 in connection with motor cable with hand control GD673!

- Use skin abrading handpiece GB280 only with motor cable GD672 and foot control.

WARNING Risk of injury and material damage due to inappropriate use of handpieces or attachments!

- Always follow the safety advice and information given in the instructions for use of the handpieces and attachments.
- Do not exceed the specified maximum motor speed.

WARNING Risk of injury caused by loose milling heads when operating high-speed motors in left rotation mode!

- Do not use screw-coupled milling heads.

CAUTION Damage to motor and handpiece!

- Couple/uncouple handpieces only with the motor inactive.

Note

For using saws GB128R, GB129, GB130R and motor cable with hand control GD673, we recommend guiding the saw with both hands for optimal handling. If single-handed saw handling is required, use motor cable GD672 with a foot control.

Connecting micro-Line handpieces to low-speed motors GD678 and GD679


Fig. 15

Legend

- A Handpiece
- B Catch
- C Centering lug
- 32 Motor

- To couple handpiece **A**:
Install handpiece **A** on centering lug **C**.
Handpiece **A** engages in its proper position.
- Pull at motor **32** to check that the handpiece is securely coupled.
- To uncouple handpiece **A**:
press catch **B** and remove handpiece **A**.

Connecting Hi-Line handpieces to high-speed motors GD676 and GD677


Fig. 16

Legend

- A Handpiece
- B Threaded ring
- C Recess
- D Lugs
- 32 Motor

	Symbol "Handpiece release"
	Symbol "Handpiece locking"

Handpiece **A** can be coupled to motor **32** in different positions.

Coupling:

- To release handpiece **A**: Turn threaded ring **B** in the direction indicated by the arrow in symbol "Handpiece release" as far as it will go.
- To couple handpiece **A**: Install handpiece **A** at the desired position on motor **32**, making certain that lugs **D** of motor **32** engage in recesses **C** at handpiece **A**.
- Lock handpiece **A**: Turn threaded ring **B** in the direction indicated by the arrow in symbol "Handpiece locking" as far as it will go (to prevent loosening of the handpiece due to vibrations during operation).

Uncoupling:

- To uncouple the handpiece: Turn threaded ring **B** in the direction indicated by the arrow in symbol "Handpiece release" as far as it will go and detach handpiece **A**.

microspeed® uni motor system

Connecting Hi-Line XS handpieces to high-speed motors GD676 and GD677


Fig. 17

Legend

- A Handpiece
- B Threaded ring
- C Handpiece coupling
- 32 Motor

	Symbol "Handpiece release"
	Symbol "Handpiece locking"

Coupling:

- Turn threaded ring **B** of handpiece **A** in the direction indicated by the arrow in symbol "Handpiece release" as far as it will go.
- Install handpiece **A** on handpiece coupling **C** of motor **32** in such a way that the pin on the handpiece engages in one of the notches of the handpiece coupling.
- Turn threaded ring **B** in the direction indicated by the arrow in symbol "Handpiece locking" as far as it will go.
- Make certain that threaded ring **B** is tightened to its limit stop to prevent it from coming loose due to vibrations during operation.

Uncoupling:

- Turn threaded ring **B** in the direction indicated by the arrow in symbol "Handpiece release" as far as it will go.
- Detach handpiece **A** from motor **32**.

Connecting mini-Line attachments to microspeed® uni mini 100 motor GD674 and microspeed® uni mini pistol handpiece GD684


Fig. 18

Legend

- A mini-Line attachment
- B Tension ring
- C Motor drive

Coupling:

- Turn tension ring **B** on motor **C** in the "open" direction as indicated by the arrow, as far as it will go.
- Push the coupling of mini-Line attachment **A** into tension ring **B** at motor **C** as far as it will go.
- Turn tension ring **B** on motor **C** against the "open" direction as indicated by the arrow, as far as it will go.

Uncoupling:

- Turn tension ring **B** on motor **C** in the "open" direction as indicated by the arrow, as far as it will go.
- Retract mini-Line attachment out of tension ring **B** at motor **C**.
- If tension ring **B** on the coupling cannot be loosened by hand, use the auxiliary wrench of motor **C**.

Connecting protective sleeve GD684225 to microspeed® uni mini pistol handpiece GD684

Fig. 19

Legend

A Protective sleeve GD684225

33 Pistol handpiece

Coupling:

- Slide protective sleeve **A** onto pistol handpiece **33** so that it clicks into position.

Uncoupling:

- Retract protective sleeve **A** from pistol handpiece **33**.

Coupling the blades to shaver handpiece GD682

Risk of injury and material damage due to inappropriate use of the blades!

- WARNING** ➤ Always follow the safety advice and information given in the respective instructions for use.

Note

Only use Aesculap blades.

Coupling a blade to the shaver handpiece


Fig. 20

Legend

A Blade

B Cap

C Fixation nose

D Coupling

E Unlocking sleeve

F Groove

37 Shaver handpiece

Blade **A** can be mounted on shaver handpiece **37** in four different positions.

- Hold blade **A** at cap **B** and push it into the blade adapter of the handpiece as far as it will go, making certain that fixation nose **C** sits in groove **F** of the motor.

The blade window is aligned with fixation nose **C**.

- Check the firm seating of blade **A** by pulling blade **A**.

microspeed® uni motor system

Uncoupling a blade from the shaver handpiece


Fig. 21
 ➤ Pull back unlocking sleeve E at the handpiece in the direction of the arrow; remove blade A, holding it at cap B.

Connecting a suction tube to shaver handpiece GD682


Fig. 22

Legend

- A Suction tube
- 36 Suction control lever
- 37 Shaver handpiece
- 38 Tube olive

- To connect suction tube A to shaver handpiece 37: Push suction tube A (GF572 or GF573) of the surgical suction device over tube olive 38.

4.2 Function checks

- Prior to each use, inspect the entire unit for functionality and proper condition.
- Verify that all components to be used function properly and are properly connected.

The connected motor type must be indicated in the appropriate motor control field on the display.

- Set the motor release switch to its "ON" position, if necessary.
 The motor to be used is released.
- Actuate the foot or hand control down to the limit stop.
 The motor starts and reaches the maximum speed indicated in the motor control field in display 2.
 The motor audibly runs smoothly and at a constant speed.
 The qualitative display of the actual motor speed in the motor control field is fully illuminated.
- Enable coolant pump 4 in the respective motor control field or via the function key of the foot or hand control device.

Coolant pump 4 is started by activating the motor.

Once the motor was active for more than two seconds, a special drive control prevents further coolant from dripping out.

4.3 Safe operation

Risk of injury caused by inadvertent activation of the motor!

- Secure any motors that are not actively operated against inadvertent activation through foot or hand control.

Risk of burns to patient and user, due to hot motor/hot handpiece/hot tool!

- Put down the hot motor/hot handpiece/hot tool beyond reach of the patient.

Risk of injury and material damage due to inappropriate use of handpieces or motors!

- Always follow the safety advice and information given in the instructions for use of the handpieces and motors.
- Do not exceed the specified maximum motor speed.

Risk of injury caused by loose milling heads when operating high-speed motors in left rotation mode!

- Do not use screw-coupled milling heads.

4.3.1 Setting up and operating the motors GD674, GD676, GD677, GD678 and GD679

Note

The instruction menu allows changing motor type-specific settings (run-down and run-up ramp), see chapter Instruction menu.

A motor can be operated and settings can be changed only if:

- a motor **32** is connected at control unit **1**
- the motor control field is not locked, meaning
 - when using the motor cable for foot control GD672, the motor release switch is in its "ON" position
 - when using the motor cable for foot control GD672, a foot control device is connected at control unit **1**
- any motor connected to the other motor connection socket is inactive.

The following settings cannot be changed while motor **32** is running:

- Motor direction (of rotation)
- Motor speed range upper limit
- Operating mode

All settings for the respective motor type and motor connection socket are saved on exit from the settings mode.

Changing the settings in the Device settings menu on the control unit

Fig. 23

Selecting the motor speed range upper limit

The set motor speed upper limit is indicated in the respective motor control field in display **2** and can be adjusted via the touch control keys.

Note

The setting range and the step width for the speed range upper limit depend on which motor was selected, see chapter Technical specifications.

- To increase the motor speed range upper limit: press key to increase the motor speed range upper limit.
- To reduce the motor speed range upper limit: press key to reduce the motor speed range upper limit.

Selecting the saw program (only for motors GD678 and GD679)

This key can only be activated when no motor is active.

The saw program is activated and deactivated via the control key in the respective motor control field.

The same control key also serves as a status indicator:

- The symbol displayed on a white background means that the saw program is not active.
- The symbol displayed on a green background means that the saw program is active, with a maximum setting of 16,000 1/min for the motor speed range upper limit.

The settings most recently used with the saw program (motor speed range upper limit, dose level for coolant pump **4**) are recalled automatically.

microspeed® uni motor system

Setting the coolant pump delivery rate

Coolant pump 4 is activated and deactivated via the control key in the respective motor control field.

The same control key also serves as a status indicator:

- The symbol displayed on a white background means that the pump is switched off.
- The symbol displayed on a green background means that the pump is switched on.
- To increase the delivery rate: press the key to increase the delivery rate in the pump control field.
- To reduce the delivery rate: press the key to reduce the delivery rate in the pump control field.

Changing the settings in the Device settings menu through the foot control

- Press and hold down foot control function button 18 for three seconds. Control unit 1 switches to the settings mode. Motor 32 and coolant pump 4 are disabled.

Saw program selection (only for motors GD678 and GD679)

Notes

These motors support an additional saw program. The maximum motor speed and the coolant delivery rate are changed in the same way as for motor types GD674, GD676 and GD677.

In the saw program, the maximum motor speed is limited to 16,000 1/min. The settings that were used most recently with or without the saw program are recalled and displayed automatically.

The indicator/control key for the saw program flashes in the respective motor control field.

The saw program is activated or deactivated by pressing the respective pedal of foot control 21/24:

- Symbol on white background: program deactivated
- Symbol on green background: program activated
- To switch from parameter saw program to parameter motor speed range upper limit: briefly press foot control function button 18.

Changing the motor speed range upper limit

The motor speed range upper limit display for the motor is flashing in the motor control field.

- To increase the motor speed range upper limit:
 - press right pedal 19 when using foot control GD671.
 - for foot control GD668, set motor direction switch 22 to right rotation and press the pedal.

If the motor speed range upper limit is reached and the pedal is still held down, the display switches to the minimum setting for the motor speed range upper limit before rising according to the technical specifications of the respective motor type.

- To reduce the motor speed range upper limit:
 - press left pedal 20 when using foot control GD671.
 - for foot control GD668, set motor direction switch 22 to left rotation and press the pedal.

If the minimum possible setting for the motor speed range upper limit is reached and the pedal is still held down, the display switches to the maximum setting for the motor speed range upper limit before being reduced according to the technical specifications of the respective motor type.

- To switch from parameter motor speed range upper limit to parameter coolant pump 4: briefly press foot control function button 18.

Setting the coolant pump delivery rate

The display for the selected dosage is flashing in the pump control field.

- To increase the delivery rate:
 - press right pedal 19 when using foot control GD671.
 - for foot control GD668, set motor direction switch 22 to right rotation and press the pedal.

If the maximum dosage is reached while the pedal is still held down, the display switches to the minimum dosage setting and then increases according to the preselected dosage.

- To reduce the delivery rate:
 - press left pedal 20 when using foot control GD671.
 - for foot control GD668, set motor direction switch 22 to left rotation and press the pedal.

If the minimum dosage is reached and the pedal still held down, the display switches to the maximum dosage setting before falling back to the preselected dosage.

Exiting the settings mode

- Press and hold down foot control function button 18 for three seconds. The settings mode is quitted automatically if no button or pedal is pressed within five seconds.

Changing the settings in the Device settings menu through hand control

- Press and hold down hand control function button 26 for three seconds.

Control unit 1 switches to the settings mode.

Motor 32 and coolant pump 4 are disabled.

Selecting the saw program (only for motors GD678 and GD679)

Notes

These motors support an additional saw program. The maximum motor speed and the coolant delivery rate are changed in the same way as for motor types GD674, GD676 and GD677.

In the saw program, the maximum motor speed is limited to 16,000 1/min. The settings that were used most recently with or without the saw program are recalled and displayed automatically.

The indicator/control key for the saw program flashes in the respective motor control field.

The saw program is activated and deactivated through lever **29**:

- Symbol on white background: program deactivated
- Symbol on green background: program activated
- To switch from parameter saw program to parameter motor speed range upper limit: briefly press hand control function button **26**.

Changing the motor speed range upper limit

The motor speed range upper limit display for the motor is flashing in the motor control field.

- To increase/reduce the motor speed range upper limit: actuate lever **29**. If the motor speed upper limit is reached while lever **29** is still being actuated, the display switches to the minimum setting for the motor speed upper limit before rising according to the technical specifications of the respective motor type.
- To switch from parameter motor speed range upper limit to parameter coolant pump: briefly press hand control function button **26**.

Setting the coolant pump delivery rate

The display for the selected dosage is flashing in the pump control field.

- To increase or decrease the delivery rate: Actuate the lever until the required delivery rate is displayed. As soon as the maximum delivery rate is reached, the lowest rate will be displayed again.
- To switch from parameter coolant pump to parameter motor direction: briefly press hand control function button **26**.

Setting the motor direction

The most recently set motor direction flashes in the motor control field.

- Actuate lever **29** until the indicator for the required motor direction lights up.

Exiting the settings mode

- Press and hold down hand control function button **26** for three seconds. The settings mode is automatically exited if no key or control element is actuated within five seconds.

Activating a motor through foot control

To activate right rotation:

- press right pedal **19** when using foot control GD671.
- when using foot control GD668, set motor direction switch **22** to right rotation and press pedal **23**.

The motor direction indicator for right rotation is illuminated in the motor control field and motor **32** runs clockwise.

To activate left rotation:

- press left pedal **20** when using foot control GD671.
- For foot control GD668, set motor direction switch **22** to left rotation and press pedal **23**.

The motor direction indicator for left rotation is illuminated in the motor control field and motor **32** runs counterclockwise.

Control unit **1** issues an acoustic signal.

- To switch on/off coolant pump **4**: briefly press foot control function button **18**.

Activating a motor with a handpiece

The selected motor direction is displayed in the motor control field.

- Actuate lever **29**. Motor **32** runs in the direction indicated in the motor control field.
- To switch on/off coolant pump **4**: briefly press hand control function button **26**.

microspeed® uni motor system

4.3.2 Setting up and operating the microspeed® uni shaver handpiece GD682

Risk of injury caused by inadvertent activation of the motor!

- Secure any motors that are not actively operated against inadvertent activation through foot or hand control.

Risk of injury and material damage due to inappropriate use of the blades!

- Always follow the safety advice and information given in instructions for use TA011579.
- Do not exceed the specified maximum motor speed.

Note

Shaver handpiece GD682 can only be used with foot control GD671 and a surgical suction device (e.g. GF200 or GF210).

The instruction menu allows changing motor type-specific settings (slow-down and speed-up ramps and oscillation frequency, speed-up ramp for Shaving mode), see chapter Instruction menu.

Parameter settings can only be changed if:

- a motor with motor cable GD673 is connected at the second motor connection socket without being active, or
- a pistol handpiece GD684 is connected at the second motor connection socket without being active, or
- a motor with motor cable GD672 is connected at the second motor connection socket and the motor release switch is in its "OFF" position, or
- only one shaver handpiece GD682 is connected at the control unit and
- a foot control device is connected to the control unit.

The following settings cannot be changed while shaver handpiece GD682 is running:

- Motor direction (of rotation)
- Motor speed range upper limit
- Operating mode

All settings will be saved on exit from the settings mode.

Operating modes

There are two types of blades for shaver handpiece GD682:

- Reamer blades
- Shaver blades

Control unit 1 recognizes whether the coupled blade is a reamer or a shaver blade. The respective mode of operation (reaming or shaving) is indicated on the display.

An acoustic signal is issued whenever a blade of a different type is inserted, and the respective symbol is displayed in the appropriate motor control field.

With a reamer blade coupled, the motor can be operated in right or left rotation, but not in oscillating mode. The maximum speed is limited to 8,000 1/min. The symbol for the "Reaming" operating mode is displayed in the control field.

With a shaver blade coupled, the motor can be operated in right or left rotation as well as in oscillating mode. The maximum speed is limited to 5,000 1/min. The symbol for the "Shaving" operating mode is displayed in the control field.

Changing the settings in the Device settings menu on the control unit

Notes

For shaver handpiece GD682, the tool speed is indicated in the display.

The setting range and the step width for the speed range upper limit depend on which operating mode (cutting or shaving) is selected, see chapter Technical specifications.

Selecting the motor speed range upper limit

The set motor speed range upper limit is indicated in the respective motor control field in display 2 and can be adjusted via the touch control keys.

- To increase the motor speed range upper limit: press key to increase the motor speed range upper limit.
- To reduce the motor speed range upper limit: press key to reduce the motor speed range upper limit.

Setting the coolant pump delivery rate

Note

The coolant pump serves exclusively for cooling or irrigating the tools. The delivery rate at the highest dose level is approx. 80 ml per minute.

Coolant pump 4 is activated and deactivated via the control key in the respective motor control field.

The same control key also serves as a status indicator:

- The symbol displayed on a white background means that the pump is switched off.
- The symbol displayed on a green background means that the pump is switched on.
- To increase the delivery rate: press the key to increase the delivery rate in the pump control field.
- To reduce the delivery rate: press the key to reduce the delivery rate in the pump control field.

Changing the settings in the Device settings menu through foot control GD671

- Press and hold down foot control function button 18 for three seconds. Control unit 1 switches to the settings mode. Shaver handpiece 37 and coolant pump 4 are disabled.

Changing the motor speed range upper limit

The motor speed range upper limit display for the motor is flashing in the motor control field.

- To increase the motor speed range upper limit: press right pedal **19**.
If the motor speed range upper limit is reached and the pedal is still held down, the display switches to the minimum setting for the motor speed range upper limit before rising according to the technical specifications of the respective motor type.
- To reduce the motor speed range upper limit: press left pedal **20**.
If the minimum possible setting for the motor speed range upper limit is reached and the pedal is still held down, the display switches to the maximum setting for the motor speed range upper limit before being reduced according to the technical specifications of the respective motor type.
- To switch from parameter motor speed range upper limit to parameter coolant pump: briefly press function button foot control **18**.

Setting the coolant pump delivery rate

The selected dosage is shown as a flashing display in the pump control field.

- To increase the delivery rate: press right pedal **19**.
If the maximum dosage is reached while the pedal is still held down, the display switches to the minimum dosage setting and then increases according to the preselected dosage.
- To reduce the delivery rate: press left pedal **20**.
If the minimum dosage is reached and the pedal still held down, the display switches to the maximum dosage setting before falling back to the preselected dosage.

Exiting the settings mode

- Press and hold down foot control function button **18** for three seconds.
The settings mode is automatically exited if no key or control element is actuated within five seconds.

Activating a shaver handpiece through foot control

Notes

With the "Reaming" operating mode selected, the handpiece can be operated in right or left rotation, but not in oscillating operation.

- To activate right rotation: press the right pedal.
The motor direction indicator for right rotation is illuminated in the motor control field and shaver handpiece **37** runs clockwise.
- To activate left rotation: press the left pedal.
The motor direction indicator for left rotation is illuminated in the motor control field and shaver handpiece **37** runs counterclockwise.
Control unit **1** issues an acoustic signal.
- To switch on/off coolant pump **4**: briefly press foot control function button **18**.

Note

With the "Shaving" operating mode selected, the handpiece can be operated in right and left rotation as well as in oscillating operation.

- To activate right rotation: press the right pedal.
The motor direction indicator for right rotation is illuminated in the motor control field and shaver handpiece **37** runs clockwise.
- To activate left rotation: press the left pedal.
The motor direction indicator for left rotation is illuminated in the motor control field and shaver handpiece **37** runs counterclockwise.
- To activate oscillating operation: press both pedals at the same time.
The motor direction indicators for left and right rotation in the motor control field light up alternately.
Shaver handpiece **37** alternately runs clockwise and counterclockwise.
- To switch on/off coolant pump **4**: briefly press foot control function button **18**.

microspeed® uni motor system

Setting the suction control device

Fig. 24

Legend

36 Suction control lever

The suction power control is continuously adjustable through the suction control lever.

In position "0", suction is interrupted between the blade and the suction device (no fluid drainage via the blade).

In position "I", suction is effective at full power according to the setting.

Note

The blade can get blocked at low suction power. This can be remedied by setting the device to full suction power for a short period.

➤ Adjust suction control lever 36 to the required suction power.

4.3.3 Setting up and operating the microspeed® uni mini pistol handpiece GD684

Risk of injury and damage due to long Kirschner wires!

➤ Use a sleeve when working with long Kirschner wires.

Notes

The pistol handpiece can be locked against inadvertent actuation by gently pressing and turning the upper trigger by 90°.

The pistol handpiece cannot be used with a foot control.

Working with an attached saw handpiece is only possible in the standard operating mode.

The instruction menu allows changing motor type-specific settings (slow-down and speed-up ramps, selection of oscillating or tapping mode), see chapter Instruction menu.

The motor can only be operated if:

- any motor connected to the other motor connection socket is inactive.

The following settings cannot be changed while motor 32 is running:

- Motor speed range upper limit
- Operating mode

All settings will be saved on exit from the settings mode.

Pistol handpiece GD684 supports two operating modes:

- Standard operating mode
- Oscillating or tapping mode

Changing the settings in the Device settings menu on the control unit

Note

The setting range and the step width for the speed range upper limit depend on the selected operating mode, see chapter Technical specifications.

The speed range upper limit can be changed only in the standard operating mode and only through the control unit.

The coolant pump can be activated, and the delivery rate can only be changed, via the control unit.

Selecting the motor speed range upper limit

Note

The motor speed range upper limit can only be selected in the standard operating mode.

The set motor speed range upper limit is indicated in the respective motor control field in display 2 and can be adjusted via the touch control keys.

- To increase the motor speed range upper limit: press key to increase the motor speed range upper limit.
- To reduce the motor speed range upper limit: press key to reduce the motor speed range upper limit.

Setting/activating the coolant pump delivery rate

Coolant pump 4 is activated and deactivated via the control key in the respective motor control field.

The same control key also serves as a status indicator:

- The symbol displayed on a white background means that the pump is switched off.
 - The symbol displayed on a green background means that the pump is switched on.
- To increase the delivery rate: press the key to increase the delivery rate in the pump control field.
- To reduce the delivery rate: press the key to reduce the delivery rate in the pump control field.

Selecting the operating mode

Standard operating mode

Pistol handpiece GD684 runs in standard mode if:

- the control unit recognizes that a pistol handpiece GD684 is connected,
- "oscillating or tapping mode" was selected in the Instruction menu, but not activated through the pistol handpiece.

The symbol "Oscillating or tapping mode" is displayed in the motor control field.

Oscillating or tapping mode

Note

In this operating mode, the motor speed range upper limit cannot be changed via the touch control field on the control unit.

On system power-up, the oscillating or tapping mode is by default deactivated and has to be activated again through the pistol handpiece.

The preselected setting in the Instruction menu is saved and remains in memory even through subsequent system restarts.

Pistol handpiece GD684 will operate in "Oscillating or tapping mode" if:

- "Oscillating or tapping mode" was selected in the Instruction menu,
- this mode is activated through the pistol handpiece.

The symbol "display for oscillating or tapping mode active" is illuminated in the motor control field.

Selecting Oscillating or tapping mode in the Instruction menu

- Call up Instruction menu in display 2: press control key "Instruction menu" (arrow A).


Fig. 25

Legend

- A Arrow


Fig. 26

Legend

- A Arrow

microspeed® uni motor system

- Call up menu "Motor settings": press control key "Motor" (arrow A) in the "Settings" field.


Fig. 27

Legend

A Arrow

- Call up menu "Setting" for pistol handpiece GD684: press control key "GD684" (arrow A).


Fig. 28

Legend

A Display for oscillating or tapping mode

B Control keys for Oscillating or tapping mode "YES"

C Control key "Save setting"

The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Once the change has been saved by pressing this control key, the control key will be displayed on a white background.

D Control key "EXIT"

To exit from the Instruction menu

- To pre-select "Oscillating or tapping mode": press control key B.
- Save selection: press control key C.
- Exit the Instruction menu: press control key D.

- The symbol "Display for oscillating or tapping mode" (arrow A) is illuminated in the motor control field.

Fig. 29

Legend

A Arrow

Activating the pistol handpiece

Fig. 30

Legend

- 33 Pistol handpiece
- 34 Upper trigger
- 35 Lower trigger

Standard operating mode

- To activate right rotation: press lower trigger **35**.
The motor speed is proportional to the actual trigger travel.
The motor direction indicator for right rotation is illuminated in the motor control field.
Pistol handpiece **33** runs clockwise.
- To activate left rotation: First pull upper trigger **34**, then lower trigger **35**.
The motor direction indicator for left rotation is illuminated in the motor control field.
Pistol handpiece **33** runs counterclockwise.
Control unit **1** issues an acoustic signal.
- To reverse the motor direction: press upper trigger **34** with lower trigger **34** already pressed.
The motor direction indicator for right rotation is illuminated in the motor control field.
Pistol handpiece **33** runs clockwise.
The direction indicator for left rotation is illuminated and pistol handpiece **33** reverses its direction of rotation.
Control unit **1** issues an acoustic signal.

Oscillating or tapping mode

- For selecting the Oscillating or tapping mode in the Instruction menu, see chapter Selecting Oscillating or tapping mode in the Instruction menu.
- To activate the Oscillating or tapping mode: press and hold down upper trigger **34** for three seconds.
Control unit **1** issues an acoustic signal.
The symbol "Display for oscillating or tapping mode active" is illuminated in the motor control field.

Fig. 31

- The keys for increasing or reducing the speed range upper limit are faded out as, in this operating mode, the speed range upper limit cannot be changed via the touch control field on control unit **1**.
- To activate left rotation: press upper trigger **34** with lower trigger **35** already pressed.
The motor direction indicator for left rotation is illuminated in the motor control field.
Pistol handpiece **33** runs counterclockwise.
Control unit **1** issues an acoustic signal.
- To activate tapping: press lower trigger **35** within the first half of its total travel range.
The motor direction indicators for left and right rotation in the motor control field light up alternately.
Pistol handpiece **33** alternately runs counterclockwise and clockwise, with the clockwise angle of rotation always larger than the counterclockwise angle of rotation.
- To activate oscillating operation: press lower trigger **35** within the second half of its total travel range.
The motor direction indicators for left and right rotation in the motor control field light up alternately.
Pistol handpiece **33** runs counterclockwise and clockwise alternately.

microspeed® uni motor system

- To deactivate the Oscillating or tapping mode: press and hold down upper trigger 34 for three seconds. Control unit 1 issues an acoustic signal.

Fig. 32
The symbol "Display for oscillating or tapping mode" is illuminated in the motor control field.

4.4 Settings in the Instruction menu

Note
The Instruction menu can only be called up when no motor is being operated.

Motor operation will be blocked while the Instruction menu is active.

The Instruction menu is divided into two categories:

- Information:
Call-up of information on the error conditions issued on the display (error list) and on the software version installed
 - Setting:
To change specific settings; subdivided into the sections "Device settings" and "Motor" (type).
All settings can be reset to the factory settings by pressing the control key "Factory setting".
- To call up the Instruction menu:
press control key "Instruction menu".

- To call up submenus:
press the respective control keys.


Fig. 33

Legend

- A Section "Information"
- B Control key "Info error list":
Calling up the submenu "Information on the error conditions issued on the display"
- C Display of "Software version"
- D Section "Setting"
- E Control key "Device setting":
Calling up the submenu "Device-specific settings"
- F Control key "Motor":
Calling up the submenu "Motor (type)-specific settings"
- G Control key "Factory setting":
Reset, restoring factory settings (see chapter Factory settings)
- H Control key "EXIT":
To exit from the Instruction menu

4.4.1 Submenu information on the error conditions shown on the display

- Press control key "Info error list" in menu "Information/Setting".


Fig. 34

Legend

- A** Text field:
In this field the chapter "Troubleshooting list", Malfunctions with error numbers, is displayed
- B** Control key "Scroll up":
press this key to scroll up through the text.
If this key is pressed at the beginning of the text, the display jumps to the end of the text
- C** Control key "Scroll down":
press this key to scroll down through the text.
If this key is pressed at the end of the text, the display jumps to the beginning of the text
- D** Control key "Change to higher menu level"
(to submenu "Device-specific settings")
- E** Control key "EXIT":
To exit from the Instruction menu

4.4.2 Submenu Device-specific settings

Device-specific settings are:

- Language setting
- Volume of acoustic alarm/notification signal
- Press control key "Device setting" in menu "Information/Setting".


Fig. 35

Legend

- A** Control key "Language":
Call up submenu for language setting
- B** Control key "Volume":
Call submenu for setting the volume of the acoustic alarm/notification signal
- C** Control key "Change to higher menu level":
(to Information/Setting menu)
- D** Control key "EXIT":
To exit from the Instruction menu

Language setting

➤ Press control key "Language" in submenu "Setting".


Fig. 36

Legend

- A** Control key for language display:
Select language; the selected language will be displayed on a green background
- B** Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.
- C** Control key "Change to higher menu level"
(to submenu Device-specific settings)
- D** Control key "EXIT":
To exit from the Instruction menu

Setting the volume of acoustic alarm/notification signals

➤ Press control key "Setting" in submenu "Device setting".


Fig. 37

Legend

- A** Qualitative bar display of the volume to be set
- B** Control keys "Change volume"
- C** Control key "Volume":
To call up an acoustic signal at the presently set volume
- D** Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.
- E** Control key "Change to higher menu level"
(to submenu Device-specific settings)
- F** Control key "EXIT":
To exit from the Instruction menu

4.4.3 Submenu Motor (type)-specific settings

This submenu allows specific settings for every motor type described in the present instructions for use.

- Press control key "Motor" in menu "Information/Setting".


Fig. 38

Legend

- A** Control keys "Motor type":
Select a motor type
- B** Control key "Change to higher menu level":
(to Information/Setting menu)
- C** Control key "EXIT":
To exit from the Instruction menu

Motor (type)-specific settings GD674, GD676, GD677, GD678, GD679

For these motor types, the slow-down and speed-up ramps can be adjusted.

- Press control key "Motor" (e.g. GD678).


Fig. 39

Legend

- A** Symbol for "Speed up" ramp
- B** Qualitative bar display of the speed-up ramp to be set
- C** Control keys for changing the speed-up ramp
- D** Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.
- E** Symbol for "Slow down" ramp
- F** Qualitative bar display of the slow-down ramp to be set
- G** Control keys for changing the slow-down ramp
- H** Control key "EXIT":
To exit from the Instruction menu
- I** Control key "Change to higher menu level":
(to submenu Motor (type)-specific settings)

microspeed® uni motor system

Motor (type)-specific settings shaver handpiece GD682

The following settings are available for the shaver handpiece:

- Slow-down and speed-up ramps
- Oscillation frequency (= duration of left-right rotation, in seconds) and speed-up ramp in oscillating operation (for Shaving mode)

➤ Press control key for motor type GD 682.

Adjust slow-down and speed-up ramps:


Fig. 40

Legend

- A Symbol for "Speed up" ramp
- B Qualitative bar display of the speed-up ramp to be set
- C Control keys for changing the speed-up ramp
- D Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.
- E Symbol for "Slow down" ramp
- F Qualitative bar display of the slow-down ramp to be set
- G Control keys for changing the slow-down ramp
- H Control key "EXIT":
To exit from the Instruction menu
- I Control key "Change to higher menu level":
(to submenu Motor (type)-specific settings)
- J Control key "MORE":
Call up of submenu "Set oscillation frequency and speed-up ramp in oscillating operation" (for Shaving mode of operation)

Setting the oscillation frequency and speed-up ramp for oscillating operation in the "Shaving" mode:


Fig. 41

Legend

- A Symbol for oscillating operation
- B Control keys for changing the oscillation frequency/speed-up ramp
- C Oscillation frequency display
- D Symbol for "Speed up" ramp
- E Qualitative bar display of the speed-up ramp to be set
- F Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Press this key to save the changed setting. Once the change has been saved, the control key will be displayed on a white background.
- G Control key "Change to higher menu level":
(to submenu Motor(type)-specific settings/Motor GD682/Set the slow-down and speed-up ramps)
- H Control key "EXIT":
To exit from the Instruction menu

Motor (type)-specific settings pistol handpiece GD684

The following settings are available for the pistol handpiece:

- Slow-down and speed-up ramps
 - Oscillating or tapping mode "Yes/No"
 - Rotation angle in Oscillating or tapping mode
- Press control key "Motor" (motor type GD684) in menu "Information/Setting".


Fig. 42

Legend

- A Symbol for "Speed up" ramp
- B Symbol for "Slow down" ramp
- C Symbol for Oscillating or tapping mode
- D Control keys for oscillating mode, "Yes" and "No"
On activation of "Yes" the control key "MORE" will be displayed. If "Yes" is selected and saved via control key E, the symbol "Oscillating or tapping mode" will be displayed in the motor control field after exiting the Instruction menu, if a pistol handpiece is connected.
- E Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.
- F Control key "Change to higher menu level":
(to submenu Motor (type)-specific settings)
- G Control key "EXIT":
To exit from the Instruction menu
- H Control key "MORE":
Call-up of submenu "Set rotation angle in oscillating or tapping mode"
- I Qualitative bar display of the slow-down ramp to be set
- J Control keys for changing the slow-down ramp

- K Control keys for changing the speed-up ramp
 - L Qualitative bar display of the speed-up ramp to be set
- To set the rotation angle for Oscillating or tapping mode:
➤ Press control key "MORE".


Fig. 43

Legend

- A Symbol for settings for oscillating operation:
Changing the angle of rotation for oscillating operation.
- B Field for settings for tapping operation
In this field the tool rotation angle for tapping can be changed. Angle adjustment for right rotation (forward) and angle adjustment for left rotation (backward).
- C Symbol for right rotation
- D Symbol for left rotation
- E Control keys for changing the angle of rotation
- F Display of the angle of rotation to be set
- G Control key "Change to higher menu level":
(to submenu Motor(type)-specific settings/Motor GD684/Set slow-down and speed-up ramps, Oscillating or tapping mode Yes/No)
- H Control key "EXIT":
To exit from the Instruction menu
- I Control key "Save setting":
The control key "Save setting" will be highlighted by a green background whenever a setting was changed. Pressing this key will cause the change to be saved. Once the change has been saved, the control key will be displayed on a white background.

5. Processing

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net.

Note

Processing of accessories such as handpieces, tools and tubing is described in the respective instructions for use and product care advisories. Do not process products that are intended for single use only.

5.1 Processing the microspeed® uni control unit with cooling unit GD670 and microspeed® uni foot control devices GD668/GD671

Damage to, or destruction of the product caused by mechanical cleaning/disinfecting!

- Always clean and disinfect the product manually, never in a machine.
- Do not sterilize the product under any circumstances.

5.1.1 Cleaning/Disinfection

Risk of electric shock and fire!

- Unplug the device before cleaning it.
- Do not use flammable or explosive cleaning and disinfecting agents.
- Make certain that no fluids will penetrate the product.

Damage to the product due to inappropriate cleaning or disinfecting agents!

- Only use a cleaning/disinfecting agent that is approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.
- Use cleaning/disinfecting agents that are suitable for surface cleaning of the product.
- Always follow the manufacturers' instructions regarding concentration and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Never immerse the product in water or detergent.
- Never boil or mechanically disinfect the product.
- Let any fluid that penetrated the product drain out immediately.

5.1.2 Manual cleaning/disinfecting

- Apply cleaning and disinfectant solution.
- Wipe away cleaning and disinfectant solution with a sponge or cloth that has been moistened in clean water.
- Inspect surfaces, cavities, lumens and openings for visible debris. If necessary, repeat the cleaning/disinfection process.
- Dry the product with a lint-free cloth.

5.1.3 Control, care and inspection

- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged.

5.1.4 Storage

- Store reprocessed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

5.2 Processing the microspeed® uni motors GD674, GD676, GD677, GD678, GD679 and motor cables GD672 and GD673

5.2.1 Preparations

- Prior to first mechanical cleaning/disinfection: Mount ECCOS® storage device GB578R (for GD674, GD676, GD677, GD678, GD679) and GB689R (for GD672, GD673) in a suitable basket tray (e.g. JF214R).
- Mount the storage devices in such a way that the cables can be coiled in a circle with a bending radius of not less than 10 cm. Do not kink the cables.
- Separate the motor from the motor cable prior to reprocessing.
- Carry out non-fixating/NaCl-free pre-cleaning immediately after use.
- Remove all attached components such as adapters, tubes, etc.
- Reprocess the product immediately after use.
- Do not clean the product by ultrasound treatment under any circumstances.
- Put away the product in dry condition.

5.2.2 Cleaning/Disinfection

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- **Only use cleaning and disinfecting agents that do not affect softeners (e.g. silicone), according to the manufacturers' instructions for use.**
- **Do not exceed the maximum allowable cleaning temperature of 55 °C.**
- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions regarding concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Preferably use thermal disinfecting processes.

microspeed® uni motor system

5.2.3 Mechanical cleaning/disinfecting

Pre-cleaning

- Thoroughly pre-rinse under running water.

Cleaning

- Place the product in its proper position in the ECCOS® storage device immediately after use.
- Process the product in a single-chamber cleaning/disinfecting machine without ultrasound treatment. Follow the instructions provided by the manufacturer of the machine.
- Carry out the processing cycle:
 - Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 55 °C.
 - Wash the product for at least 3 min.
 - Neutralize if necessary.
 - Carry out intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: Rinse with distilled, demineralized or fully desalinated water for 5 min, at 93 °C.
 - Complete the program with a drying phase of at least 10 min at a temperature of not more than 120 °C.
- After completion of the mechanical cleaning/disinfecting cycle, inspect for visible debris all surfaces and places that are difficult to access.
- Clean manually if necessary.

5.2.4 Manual cleaning/disinfecting

- Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
- Do not clean the product by ultrasound treatment under any circumstances.
- Do not immerse the product in any fluid. Allow any fluid that has penetrated the device to run out immediately.
- Clean all surfaces with a soft nylon brush and a suitable cleaning/disinfecting agent.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Thoroughly brush and rinse surfaces under running water.
- Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect for visible debris all surfaces and places that are difficult to access. If necessary, repeat the cleaning/disinfection process.
- Dry the product with a lint-free cloth.
- Make certain that lumens and channels are dried, too.

5.2.5 Control, care and inspection

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged and does not show atypical running noises, overheating or excessive vibrations.
- Set aside the product if it is damaged.
- After cleaning and before sterilization, always lubricate the motor with oil spray STERILIT® HI GA536:
 - Fully insert the spray can adapter on the nipple at the rear end of the motor.
 - Hold the spray can in a vertical position and push down the spray nozzle valve until oil extrudes at the front end of the handpiece coupling of the motor.
 - Wipe away any excess oil.

5.2.6 Packaging

- Always follow the instructions for use of the respective storage devices (e.g. instructions for use TA009721 for the Aesculap ECCOS® storage system).
- Insert the product in its proper position in the ECCOS® storage device, or put it on a tray in such a way that the product is protected against damage. Observe the weight limit for each tray/container.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Make certain that the product is wrapped in such a way that recontamination of the product in the period between reprocessing and reuse will be prevented.

5.2.7 Sterilization method and parameters

Damage to motor cables GD672 and GD673 caused by sterilization in sterilizers of small volume! Hot wall surfaces (>150 °C) in the sterilizer can damage the motor cable!

- Protect the motor cable (e.g. by wrapping it in a cloth).

- Sterilize with steam, taking note of the following: Carry out steam sterilization via an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN ISO 17665 or EN 554/ISO 13683). In the fractioned vacuum process, sterilization must be carried out using the 134 °C/2 bar method with a minimum holding time of 5 minutes.

Note
Steam sterilization at 134 °C for a holding time of 18 min can adversely affect the service life of the product.

- Make certain that the motor cables are not kinked or coiled up too tightly. The coil radius has to be larger than 10 cm.
- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

5.2.8 Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

5.2.9 Storage

- Store processed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

5.2.10 Maintenance

To ensure reliable operation, Aesculap recommends that maintenance be performed after 300 reprocessing cycles, at least once per year. For services to this end, please contact your national B. Braun/Aesculap agency, see chapter Technical Service.

5.3 Processing the microspeed® uni shaver handpiece GD682

5.3.1 Preparations

- Prior to first mechanical cleaning/disinfection: Mount ECCOS® storage GB580R in a suitable tray (e.g. JF214R).
- Pre-clean the product, if necessary.
- Reprocess the product immediately after use.
- Do not immerse the shaver handpiece.
- Do not clean the shaver handpiece through ultrasound treatment.
- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions regarding concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/ alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Preferably use thermal disinfecting processes.

5.3.2 Mechanical cleaning/disinfecting

- Place the product in its proper position in ECCOS® storage GB580R immediately after use.
- Connect the rinsing tube of the washing machine at tube olive **38** of shaver handpiece **37**.
- Turn suction control lever **36** of shaver handpiece **37** to the I-position.
- Seal all rinsing tubes of the washing machine that are not used.
- Carry out the processing cycle:
 - Use an appropriate neutral, alkaline or mildly alkaline cleaning/ disinfecting agent according to the manufacturer's instructions.
 - Do not exceed the maximum washing temperature of 55 °C.
 - Wash the product for at least 3 min.
 - Neutralize if necessary.
 - Carry out intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: Rinse for 5 min at 93 °C, with distilled, demineralized or fully desalinated water.
 - Conclude the program with a drying period of at least 10 min at a temperature not exceeding 120 °C.
- After completion of the mechanical cleaning/disinfecting cycle, inspect the surfaces, the driver and the suction channel for visible debris.
- Clean manually if necessary.

microspeed® uni motor system

5.3.3 Manual cleaning/disinfecting

- Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions. Do not immerse the product in any fluid.
- Remove encrusted materials with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Clean difficult-to-access areas with soft round plastic brushes of fitting diameter.
- Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect surfaces and difficult-to-access areas for visible debris. If necessary, repeat the cleaning/disinfection process.
- Use a lint-free cloth or a compressed-air gun for drying the product.
- Make certain that difficult-to-access areas and blind holes are dried, too.

5.3.4 Control, care and inspection

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged.
- Check that the suction control lever can be moved.
- Spray-lubricate the friction surfaces of the unlocking sleeve and the suction control lever, using Aesculap STERILIT® oil spray.

5.3.5 Packaging

- Always follow the instructions for use of the respective packaging and storage equipment (e.g. instructions for use TA009721 for Aesculap ECCOS® storage system).
- Insert the product in its proper position in ECCOS® storage GB580R, or put it on a tray in such a way that the product is protected against damage. Observe the weight limit for each tray/container.
- Pack trays in a way that is appropriate for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

5.3.6 Sterilization method and parameters

Damage to motor cables GD672 and GD673 caused by sterilization in sterilizers of small volume! Hot wall surfaces (>150 °C) in the sterilizer can damage the motor cable!

- **Protect the motor cable (e.g. by wrapping it in a cloth).**
- Sterilize with steam, taking note of the following: Carry out steam sterilization via an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN ISO 17665 or EN 554/ISO 13683). In the fractioned vacuum process, sterilization must be carried out using the 134 °C/2 bar program with a minimum holding time of 5 minutes.

Note

Steam sterilization at 134 °C for a holding time of 18 min can adversely affect the service life of the product.

- Make certain that the motor cables are not kinked or coiled up too tightly. The coil radius has to be larger than 10 cm.
- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

5.3.7 Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

5.3.8 Storage

- Store reprocessed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

5.3.9 Maintenance

To ensure reliable operation, Aesculap recommends the shaver handpiece GD682 to be serviced after 300 processing cycles or at least once a year. For the respective services, please contact your national B. Braun/Aesculap agency, see chapter Technical Service.

5.4 Processing the microspeed® uni mini pistol handpiece GD684

5.4.1 Preparations

- Prior to first mechanical cleaning/disinfection: Mount ECCOS® storage device GB579R in a suitable tray (e.g. JF214R). Mount the storage device in such a way that the cables can be coiled in a circle with a bending radius of not less than 10 cm. Do not kink the cable.
- Carry out non-fixating/NaCl-free pre-cleaning immediately after use.
- Remove all attached components such as adapter, sleeve, etc.
- Reprocess the product immediately after use.
- Do not clean the product by ultrasound treatment under any circumstances.
- Put away the product in dry condition.

5.4.2 Cleaning/Disinfection

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- **Only use cleaning and disinfecting agents that do not affect softeners (e.g. silicone), according to the manufacturer's instructions for use.**
- **Do not exceed the maximum allowable cleaning temperature of 55 °C.**
- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions with regard to concentration, temperature and exposure time.
- Avoid fixating of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Preferably use thermal disinfecting processes.

5.4.3 Mechanical cleaning/disinfecting

Fig. 44

Legend

- A ECCOS® rinsing adapter
- B Tube
- 33 microspeed® uni mini pistol handpiece

Pre-cleaning

- For cleaning the cannulation, use special cleaning brush TA011327.

Cleaning

- Place the product in its proper position in ECCOS® storage device GB579R immediately after use.
- Process the product in a single-chamber cleaning/disinfecting machine without ultrasound treatment. Follow the instructions provided by the manufacturer of the machine.
- Engage ECCOS® rinsing adapter at the rear of microspeed® uni mini pistol handpiece **33** and connect it to the rinsing device of the single-chamber cleaning/disinfecting machine via a suitable tube.
- Carry out the processing cycle:
 - Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 55 °C.
 - Clean for at least 3 min.
 - Neutralize if necessary.
 - Apply intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: Rinse with distilled, demineralized or fully desalinated water for 5 min, at 93 °C.
 - Complete the program with a drying phase of at least 10 min at a temperature of not more than 120 °C.
- After completion of the mechanical cleaning/disinfecting cycle, inspect for visible debris all surfaces and places that are difficult to access.
- Clean manually if necessary.

microspeed® uni motor system

5.4.4 Manual cleaning/disinfecting

- Do not clean the product by ultrasound treatment under any circumstances.
- Use an appropriate neutral, alkaline or mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
- Do not immerse the product in any fluid. Allow any fluid that has penetrated the device to run out immediately.
- Clean all surfaces with a soft plastic brush and a suitable cleaning/disinfecting agent.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- For cleaning the cannulation, use special cleaning brush TA011327.
- Thoroughly brush and rinse surfaces under running water.
- Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect for visible debris all surfaces and places that are difficult to access. If necessary, repeat the cleaning/disinfection process.
- Dry the product with a lint-free cloth.
- Make certain that lumens and channels are dried, too.

5.4.5 Control, care and inspection

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged and does not show atypical running noises, overheating or excessive vibrations.
- Set aside the product if it is damaged.

5.4.6 Packaging

- Always follow the instructions for use of the respective storage devices (e.g. instructions for use TA009721 for the Aesculap ECCOS® storage system).
- Insert the product in its proper position in the ECCOS® storage device, or on a tray in such a way that the product is protected against damage. Observe the weight limit for each tray/container.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

5.4.7 Sterilization method and parameters

Damage to the pistol handpiece caused by sterilization in sterilizers of insufficient volume! Hot wall surfaces (>150 °C) in the sterilizer can damage the motor cable!

- **Protect the motor cable (e.g. by wrapping it in a cloth).**
- Sterilize with steam, taking note of the following: Carry out steam sterilization via an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN ISO 17665 or EN 554/ISO 13683). In the fractioned vacuum process, sterilization must be carried out using the 134 °C/2 bar program with a minimum holding time of 5 minutes.

Note

Steam sterilization at 134 °C for a holding time of 18 min can adversely affect the service life of the product.

- Make certain that the motor cables are not kinked or coiled up too tightly. The coil radius has to be larger than 10 cm.
- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

5.4.8 Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

5.4.9 Storage

- Store processed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

5.4.10 Maintenance

To ensure reliable operation, Aesculap recommends the product to be serviced after every 300 processing cycles or at least once a year.

For services to this end, please contact the B. Braun/Aesculap agency in your country, see chapter Technical Service.

6. Maintenance

Send in the microspeed® uni motor system with its accessories for annual inspection by the manufacturer, to the following address:

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-2700

Fax: +49 7461 16-2887

E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood, MO 63042

Aesculap Repair Hotline

Phone: +1 800 214-3392

Fax: +1 314 895-4420

7. Troubleshooting list

7.1 Malfunctions with error numbers

As soon as an error status is detected, an error window will pop up on the display and an acoustic warning will be issued.

In the event a malfunction cannot be remedied as described below:

- Have component(s) repaired by the manufacturer.

Malfunction	Message displayed	Cause	Detection/Remedy
Error 1	Error 1	Key actuation (softkey actuation): More than one key was pressed at the same time	Only press one key
Error 2	Error 2	Key (softkey) held down for too long: Control key pressed for longer than 20 s	Do not press any key for longer than 20 s
Error 3	Error 3 – Foot-Control	Pedal(s) of the foot control pressed during automatic self-test	Release pedal(s)
		Both pedals of foot control GD671 pressed	Press one pedal only Pressing both pedals is only allowed in operating mode "Shave" in conjunction with the shaver handpiece
		Plug contacts of foot control unclean	Clean the plug contacts
		Foot control defective	Replace foot control
Error 4	Error 4 – Hand-Control	Lever of motor cable GD673 or trigger of pistol handpiece GD684 actuated during automatic test	Release lever or trigger, respectively
		Motor cable GD673 defective	Replace motor cable
		Pistol handpiece GD684 defective	Replace pistol handpiece
Error 5	Error 5 – Select Motor	Speed request from foot control while the motor is locked by the motor release switch on motor cable GD672 (position "OFF")	Set motor release switch on motor cable to position "ON"
		Motor cable GD672 defective	Replace motor cable
	Error 5 – Select only one Motor	Speed request from foot control while two motors are released for operation the motor release switch on motor cable GD672 (position "ON")	Release the required motor through motor release switch (position "ON") Block the motor that is not required through motor release switch (position "OFF")
		Motor cable GD672 defective	Replace motor cable
	Error 5 – Connect only one Shaver	Speed request from foot control while two shaver handpieces GD682 are connected to the control unit	Connect only one shaver handpiece GD682 to the control unit

Malfunction	Message displayed	Cause	Detection/Remedy
Error 6	Error 6 – Motor	Speed request through foot or hand control while no motor is connected to the control unit	Connect a motor to the control unit
		Motor defective	Replace motor
		Motor cable GD672 or GD673 defective	Replace motor cable
	Error 6 – Motor/Hallsensor	Control unit detects malfunction in Hall system of the motor. Only for motors GD674, GD678, GD679, GD682 and GD684	Replace motor
		Motor cable GD672 or GD673 defective	Replace motor cable
	Error 6 – Motor blocked	Motor blocked	Stop motor actuation and release blockage
Motor defective		Replace motor	
Error 7	Error 7 – Motor overheated	Temperature: Motor has overheated, see chapter Overload cut-out	Allow motor to cool down, or replace motor
Error 8	Error 8 – Power Fail	Power supply error: The mains voltage is outside the defined range	Verify correct mains voltage
Error 9	Error 9	Current cannot be applied properly to the motor	Replace motor or motor cable, respectively
Error 10	Error 10 – Pump	Pump: Mechanical blockage of the pump	Check position of the tubing
Error 11	Error 11 – Internal Power-Supply e.g. (+35 V)	Internal supply voltage: - 35-V motor voltage monitor - +15-V voltage monitor - -15-V voltage monitor - +5-V voltage monitor	Switch on and off the control unit If the error reoccurs, have the control unit repaired by the manufacturer
Error 12	Error 12 – System Error	System error: E ² PROM, DSP, controller	
Error 13	Error 13 – Battery	Battery fault (Battery on processor card GD670801)	Have battery replaced by the manufacturer
Error 14	Error 14 – External Data-Communication	External data communication error (in connection with RS232, CAN, Ethernet)	Switch on and off the control unit If the error reoccurs, have the control unit repaired by the manufacturer
Error 15	Error 15 – Power-Circuit	Motor power stage	

microspeed® uni motor system

7.2 Other malfunctions

Malfunction	Cause	Finding	Remedy
The motor does not run	No power to control unit	Control unit not connected to mains power or not switched on (Power ON indicator dark, black display)	Connect control unit to mains power, and switch it on
	Fuses blown	Power ON indicator dark, black display	Replace fuses, see chapter Fuse change
Control unit cannot be operated The motor does not run	Internal monitoring system blocking the system	Control unit not operational although it is used correctly	Switch on and off the control unit If the error reoccurs, have the control unit repaired by the manufacturer
Connected handpiece does not run	Handpiece blocked	Motor runs when the handpiece is not connected	Have handpiece repaired by manufacturer
No coolant flow	Nozzle clogged	Coolant pump works	Clean nozzle
	Tubing inserted incorrectly	Coolant pump works	Install tubing correctly, see chapter System set-up
	Coolant pump motor defective	Coolant pump does not work	Have pump motor repaired by manufacturer
Pistol handpiece does not run	Protection against inadvertent switch-on activated	Upper trigger turned by 90°	Turn back upper trigger

7.3 Fuse change

Fatal injury from electric shock!
Open components can hold life-threatening mains voltage.
 ➤ Pull the mains plug before replacing the fuse sets.

Allowable fuses:

2 x IEC 127 - T 6.3 A breaking capacity H (1 500 A)

- Unplug the device before changing the fuses!
- Unlock lug at fuse holder **13** , using a small screwdriver.
- Remove fuse holder **13**.
- Change both fuses.
- Reinsert fuse holder **13** so that it audibly engages in its proper position.

Note

If the fuses blow frequently, the device is defective and needs to be repaired; see Technical Service.

8. Technical Service

For service, maintenance or repairs, contact your national B. Braun/ Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-2700

Fax: +49 7461 16-2887

E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood, MO 63042

Aesculap Repair Hotline

Phone: +1 800 214-3392

Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

9. Accessories/Spare parts

9.1 Accessories

9.1.1 microspeed® uni motors, motor cables and foot control devices

Designation	Art. no.
microspeed® uni mini 100 motor, mini-Line coupling	GD674
microspeed® uni Hi 150 motor, Hi-Line coupling	GD676
microspeed® uni Hi 100 motor, Hi-Line coupling	GD677
microspeed® uni micro 150 motor, micro-Line coupling	GD678
microspeed® uni micro 100 motor, micro-Line coupling	GD679
microspeed® uni shaver handpiece	GD682
microspeed® uni mini pistol handpiece, mini-Line coupling	GD684
Protective sleeve (for microspeed® uni mini pistol handpiece)	GD684225
Auxiliary wrench for mini-Line coupling	GB186209
microspeed® uni motor cable for foot control	GD672
microspeed® uni motor cable with hand control	GD673
microspeed® uni foot control, two pedals	GD671
microspeed® uni foot control, one pedal	GD668

9.1.2 Cooling unit

Designation	Art. no.
Tube set, reusable	GA385
Tube set, single-use	GA391
Tube set, single-use, for microspeed® uni	GA392SU
Spike for 500-ml infusion bottles with rubber plug	GA383R
Spike for 1 000-ml infusion bottles with rubber plug	GA382R
Spike for 1 000-ml infusion bottles with Euro thread	GA381R
Bottle holder	GD412804

9.1.3 Other components

Designation	Art. no.
microspeed® uni sterile shield for GD670	GD421
Mobile stand for microspeed® uni, adjustable height	GA411
Mobile stand for surgical devices	GD416M

microspeed® uni motor system

9.1.4 Power cord

Approved in	Color	Length (in m)	Art. no.
Europe	black	1.5	TE780
Europe	black	5.0	TE730
UK	black	5.0	TE734
USA, Canada, Japan	gray	3.5	TE735

9.1.5 Maintenance, processing and service

Designation	Art. no.
Drip feed oiler STERILIT® M	GA059
Oil spray STERILIT® Hi	GA536
Oil spray STERILIT® M	GB149
Oil spray STERILIT® mini	GB372
ECCOS® mount for microspeed® uni motors	GB578R
ECCOS® mount for microspeed® uni mini pistol handpiece	GB579R
ECCOS® mount for microspeed® uni shaver handpiece	GB580R
ECCOS® mount for microspeed® motor cable	GB689R
Cleaning brush	TA011327
Service manual microspeed® uni/arthro motor system	TA022090
Service manual microspeed® uni motors and motor cables	TA011909
Service manual microspeed® uni shaver handpiece	TA011910
Service manual microspeed® uni mini pistol handpiece	TA011912

9.2 Spare parts

Designation	Art. no.
Fuse set (melting) T 6.3 A	TA021473

10. Technical specifications

10.1 Classification acc. to directive 93/42/EEC

Art. no.	Designation	Class
GD668	microspeed® uni foot control, one pedal	I
GD670	microspeed® uni control unit with cooling unit	IIa
GD671	microspeed® uni foot control, two pedals	I
GD672	microspeed® uni motor cable for foot control	I
GD673	microspeed® uni motor cable with hand control	I
GD674	microspeed® uni mini 100 motor, mini-Line coupling	IIa
GD676	microspeed® uni Hi 150 motor, Hi-Line coupling	IIa
GD677	microspeed® uni Hi 100 motor, Hi-Line coupling	IIa
GD678	microspeed® uni micro 150 motor, micro-Line coupling	IIa
GD679	microspeed® uni micro 100 motor, micro-Line coupling	IIa
GD682	microspeed® uni shaver handpiece	IIa
GD684	microspeed® uni mini pistol handpiece, mini-Line coupling	IIa

10.2 microspeed® uni control unit with cooling unit GD670

Mains voltage ranges (Current intake)	100-120 V~ (4.2-3.2 A) 220-240 V~ (1.7-1.5 A)
Frequency	50-60 Hz
Protection class (acc. to IEC 60601-1)	I
Output voltage	max. 35 V
Output current	max. 8.5 A
Device fuse Breaking capacity	T 6.3 A H (1 500 A at 250 V/50 Hz)
Delivery rate of cooling pump (at maximum dosage setting)	approx. 80 ml/min
Weight	8.9 kg
Dimensions (L x W x H)	305 mm x 305 mm x 175 mm
Application component	Type BF
EMC	IEC 60601-1-2
Norms	IEC 60601-1

10.2.1 Factory settings

Parameter	Setting/Status
Motor speed range upper limit for high-speed motors GD676 and GD677	60,000 1/min
Motor speed range upper limit for low-speed motors GD678 and GD679 with micro-Line coupling	20,000 1/min
Motor speed range upper limit for low-speed motors GD678 and GD679 with micro-Line coupling	Saw program 16,000 1/min
Motor speed range upper limit for low-speed motors with mini-Line coupling GD674 and pistol handpiece GD684	16,000 1/min
Oscillating or tapping mode (for pistol handpiece GD684)	deactivated (Instruction menu)
Rotation angle, oscillating mode (for pistol handpiece GD684)	210° clockwise, then 210° counterclockwise at gear reduction 16:1
Rotation angle for tapping mode (for pistol handpiece GD684)	270° clockwise, then 90° counterclockwise at gear reduction 16:1
Motor speed range upper limit for shaver handpiece GD682, "Shaving" mode	800 1/min
Motor speed range upper limit for shaver handpiece GD682, "Reaming" mode	6,000 1/min
Oscillation frequency in "Shaving" mode	1 s left, 1 s right

Parameter	Setting/Status
Motor type-specific standard settings for speed-up ramps	Standard
Motor type-specific standard settings for slow-down ramps	Standard
Pump operation	OFF
Pump flow	2/3-dosage
Language	English
Acoustic alarm/notification signal	max. volume

10.3 Low-speed motors

10.3.1 Low-speed motors with micro-Line coupling

Art. no.	Power output (in W)	Design
GD678	150	Standard
GD679	100	Short

Motor speeds/setting

Note

The speed figures refer to motor speeds.

Parameter	Speed/Setting
Step width for speed range upper limit setting	1,000 1/min
Motor direction settings (optical and acoustic notification for left rotation)	- right rotation - left rotation
Maximum speed range upper limit - for motor GD678 - for motor GD679	- 30,000 1/min - max. 40,000 1/min
Minimum speed range upper limit	1,000 1/min
Minimum speed setting	500 1/min
Factory setting	20,000 1/min
Saw program (with limited motor speed)	max. 16,000 1/min

microspeed® uni motor system

10.3.2 Low-speed motor with mini-Line coupling

Art. no.	Power output (in W)	Design
GD674	100	Standard

Motor speeds/setting

Note
The speed figures refer to motor speeds.

Parameter	Speed/Setting
Step width for speed range upper limit setting	1,000 1/min
Motor direction settings (optical and acoustic notification for left rotation)	- right rotation - left rotation
Maximum speed range upper limit	16,000 1/min
Minimum speed range upper limit	1,000 1/min
Minimum speed setting	500 1/min
Factory setting	16,000 1/min

10.4 High-speed motors

10.4.1 High-speed motors with Hi-Line coupling

Art. no.	Power output (in W)	Design
GD676	150	Standard
GD677	100	Short

Motor speeds/setting

Note
The speed figures refer to motor speeds.

Parameter	Speed/Setting
Step width for speed range upper limit setting	5,000 1/min
Motor direction settings (optical and acoustic notification for left rotation)	- right rotation - left rotation
Maximum speed range upper limit	80,000 1/min
Minimum speed range upper limit	10,000 1/min
Minimum speed setting	5,000 1/min
Factory setting	60,000 1/min

10.5 microspeed® uni shaver handpiece GD682

Note
The speed figures refer to tool speeds.

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Note

The step width set for "to the right" (forward) must always be larger than the step width set "to the left" (backward).

Step width/Settings oscillating operation

Parameter	Setting
Step width settings	15°-angle at the tip (16:1 gear)
Setting range	45–405°
Factory setting	210° (right and left)

Note

The step width settings for rotation "to the right" (forward) and "to the left" (backward) are always identical.

11. Disposal

Note

The user institution is obliged to process the product before its disposal, see chapter Processing.

Always adhere to national regulations when disposing of or recycling the product or its components!

Products carrying this symbol are subject to separate collection of electric and electronic devices. Within the European Union, disposal is carried out without fees by the manufacturer.

If you have any questions concerning the disposal of the product, please contact your national B. Braun/Aesculap agency, see chapter Technical Service.

12. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG & Co. KG